

Special Complex Enzyme for Animal Protein Hydrolysis

Cat: C8802

Size: 500g

Storage: Store at room temperature.

Product Parameters:

Appearance(properties): white to off-white powder

solubility: 10mg/mL water(turbidity solution)

Enzyme activity: 1500U/g

Introduction:

Used in food and feed, cosmetics, nutrition and health care and other protein hydrolysis treatment industry, according to the production needs of appropriate use.

Protocols:

According to the production needs to determine the amount of enzyme, add to the substrate solution and stir, substrate concentration 25%-50%, add amount 0.1~0.3%(dry weight of the material), hydrolysis time 3~5h(depending on the degree of hydrolysis) under the reaction conditions of temperature 50°C and natural PH value.

Notes:

1. The biochemical reagents produced by our company are basically non-aseptic packaging without special marks. If used in cell experiments, please pre-treat them in advance.
2. Once it is prepared into a solution, please pack it separately and store it to avoid product failure caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. This product is for scientific research only. Do not use for medicine, clinical diagnosis or treatment, food and cosmetics. Do not store in ordinary residential areas.
5. For your safety and health, please wear a lab coat and disposable gloves and a mask.