

1M imidazole solution

Cat No.: I8092

Package: 100mL

Storage: RT

Period of validity: one year

Description of products

Imidazole is an organic compound, a diazole, a five-membered aromatic heterocyclic compound containing two internitrogen atoms in its molecular structure. Imidazole is acidic and basic, and can form salts with strong bases.

This product is 1M imidazole solution, can be used for anti-metabolism, anti-histamine; For the purification of histidine tag protein.

Product Parameter

English name: Imidazole solution

Concentration: 1M

CAS: 288-32-4

Molecular formula: C₃H₄N₂

Molecular weight: 68.08

Appearance (character) : colorless liquid

Note

- 1、 This product is mainly composed of imidazole and deionized water.
- 2、 This product is for scientific use only. Do not use for medicine, clinical diagnosis or treatment, food and cosmetics. Do not store in ordinary residential areas.
- 3、 For your safety and health, please wear a lab coat, disposable gloves and a mask.

Related products

I8090 imidazole

Note: For more information about this product, please refer to Solarbio website.