

8M Urea solution

Cat: U8022

Package: 100mL/500mL

Storage: RT, Valid for 1 year.

Product description:

Urea also known as urea, carbamide, is one of the simplest organic compounds. For protein denaturation, as a mild solvent for insoluble or denatured proteins; Together with 6M guanidine hydrochloride and DTT, the protein in the sample was denatured.

This product is 8M urea solution, which can be used directly for experiments.

Product parameters:

1. English name: urea solution (尿素溶液)
2. Concentration: 8M
3. CAS: 57-13-6
4. Molecular formula: $\text{CH}_4\text{N}_2\text{O}$
5. Molecular weight: 60.04
6. Appearance (character) : colorless liquid

Notes:

1. This product is for scientific research only. Do not use for medicine, clinical diagnosis or treatment, food and cosmetics. Do not store in ordinary residential areas.
2. For your safety and health, please wear a lab coat, disposable gloves and a mask.

Related products:

U8020 Urea

SU8070 Urea standard